

# InVu® Ceramic Brackets with Readi-Base® eXact™ Pre-Applied Adhesive

Roth Prescription

## PATIENT KITS — MINI

	Cuspid Hooks	0.46 mm (.018") Slot	0.56 mm (.022") Slot
<b>Kits contain all ceramic brackets</b>			
Upper 5-5	Standard	<b>297-060RBX</b>	<b>297-065RBX</b>
Upper 5-5	No	<b>297-107RBX</b>	<b>297-111RBX</b>
Upper & Lower 5-5	Standard	<b>297-070RBX</b>	<b>297-075RBX</b>
Upper & Lower 5-5	No	<b>297-109RBX</b>	<b>297-113RBX</b>



## MAXILLARY BRACKETS — MINI

	Slot Size	Left/Right	Bondable
Central	0.46 mm	L	<b>297-138RBX</b>
Torque +14°	(.018")	R	<b>297-139RBX</b>
Tip +5°	0.56 mm	L	<b>297-150RBX</b>
	(.022")	R	<b>297-151RBX</b>
Lateral	0.46 mm	L	<b>297-148RBX</b>
Torque +7°	(.018")	R	<b>297-149RBX</b>
Tip +8°	0.56 mm	L	<b>297-152RBX</b>
	(.022")	R	<b>297-153RBX</b>
Cuspid	0.46 mm	L	<b>297-158RBX</b>
Torque +0°	(.018")	R	<b>297-159RBX</b>
Tip +8°	0.56 mm	L	<b>297-154RBX</b>
	(.022")	R	<b>297-155RBX</b>

Package of 5

## MANDIBULAR BRACKETS — MINI

	Slot Size	Left/Right	Bondable
Anteriors	0.46 mm	L/R	<b>297-169RBX</b>
Torque -1°	(.018")		
Tip 0°	0.56 mm	L/R	<b>297-160RBX</b>
	(.022")		
Cuspid	0.46 mm	L	<b>297-176RBX</b>
Torque -11°	(.018")	R	<b>297-177RBX</b>
Tip +6°	0.56 mm	L	<b>297-174RBX</b>
	(.022")	R	<b>297-175RBX</b>

Package of 5

## PATIENT KITS

	Cuspid Hooks	0.46 mm (.018") Slot	0.56 mm (.022") Slot
<b>Kits contain all ceramic brackets</b>			
Upper & Lower 5-5	Yes	<b>297-136RBX</b>	<b>297-137RBX</b>
Upper 5-5	Yes	<b>297-124RBX</b>	<b>297-125RBX</b>

## MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46 mm	L	<b>297-140RBX</b>
Torque +14°	(.018")	R	<b>297-141RBX</b>
Tip +5°	0.56 mm	L	<b>297-142RBX</b>
	(.022")	R	<b>297-143RBX</b>
Lateral	0.46 mm	L	<b>297-144RBX</b>
Torque +7°	(.018")	R	<b>297-145RBX</b>
Tip +8°	0.56 mm	L	<b>297-146RBX</b>
	(.022")	R	<b>297-147RBX</b>
Cuspid with Distal Hook	0.46 mm	L	<b>297-222RBX</b>
Torque 0°	(.018")	R	<b>297-223RBX</b>
Tip +8°	0.56 mm	L	<b>297-224RBX</b>
	(.022")	R	<b>297-225RBX</b>
Universal Bicuspid without Hook	0.46 mm	L/R	<b>297-226RBX</b>
Torque -7°	(.018")		
Tip 0°	0.56 mm	L/R	<b>297-228RBX</b>
	(.022")		

Package of 5

## MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anteriors	0.46 mm	L/R	<b>297-170RBX</b>
Torque -1°	(.018")		
Tip 0°	0.56 mm	L/R	<b>297-172RBX</b>
	(.022")		
Cuspid with Distal Hook	0.46 mm	L	<b>297-178RBX</b>
Torque -11°	(.018")	R	<b>297-179RBX</b>
Tip +6°	0.56 mm	L	<b>297-180RBX</b>
	(.022")	R	<b>297-181RBX</b>
1st Bicuspid without Hook	0.46 mm	L/R	<b>297-182RBX</b>
Torque -17°	(.018")		
Tip 0°	0.56 mm	L/R	<b>297-184RBX</b>
	(.022")		
2nd Bicuspid without Hook	0.46 mm	L/R	<b>297-186RBX</b>
Torque -22°	(.018")		
Tip 0°	0.56 mm	L/R	<b>297-188RBX</b>
	(.022")		

Package of 5

# InVu® Ceramic Brackets with Readi-Base® eXact™ Pre-Applied Adhesive

MBT (McLaughlin, Bennet, Trevisi System) Prescription

## PATIENT KITS

	Cuspid Hooks	0.46 mm (.018") Slot	0.56 mm (.022") Slot
Upper & Lower 5-5	Yes	297-206RBX	297-208RBX
	No	297-207RBX	297-209RBX
Upper 5-5	Yes	297-194RBX	297-196RBX

Kits contain maxillary first and second bicusps with .027" in/out compensation.



## MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46 mm (.018")	L	297-210RBX
Torque +17°		R	297-211RBX
Tip +4°	0.56 mm (.022")	L	297-212RBX
		R	297-213RBX
Lateral	0.46 mm (.018")	L	297-214RBX
Torque +10°		R	297-215RBX
Tip +8°	0.56 mm (.022")	L	297-216RBX
		R	297-217RBX
Cuspid without Hook	0.46 mm (.018")	L	297-218RBX
Torque -7°		R	297-219RBX
Tip +8°	0.56 mm (.022")	L	297-220RBX
		R	297-221RBX
Cuspid with Distal Hook	0.46 mm (.018")	L	297-222RBX
Torque 0°		R	297-223RBX
Tip +8°	0.56 mm (.022")	L	297-224RBX
		R	297-225RBX
Universal Bicuspid without Hook	0.46 mm (.018")	L/R	297-226RBX
Torque -7°			
Tip 0°	0.56 mm (.022")	L/R	297-228RBX
In/Out Comp. .027"			

Package of 5

## MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anteriors Narrow	0.46 mm (.018")	L/R	297-254RBX
Torque -6°			
Tip 0°	0.56 mm (.022")	L/R	297-255RBX
Anteriors Wide	0.46 mm (.018")	L/R	297-236RBX
Torque -6°			
Tip 0°	0.56 mm (.022")	L/R	297-237RBX
Cuspid without Hook	0.46 mm (.018")	L	297-256RBX
Torque -6°		R	297-257RBX
Tip +3°	0.56 mm (.022")	L	297-258RBX
		R	297-259RBX
Cuspid with Distal Hook	0.46 mm (.018")	L	297-260RBX
Torque 0°		R	297-261RBX
Tip +3°	0.56 mm (.022")	L	297-262RBX
		R	297-263RBX
1st Bicuspid without Hook	0.46 mm (.018")	L/R	297-264RBX
Torque -12°			
Tip 0°	0.56 mm (.022")	L/R	297-266RBX
2nd Bicuspid without Hook	0.46 mm (.018")	L/R	297-272RBX
Torque -17°			
Tip 0°	0.56 mm (.022")	L/R	297-274RBX
Universal Bicuspid without Hook	0.46 mm (.018")	L/R	297-234RBX
Torque -7°			
Tip 0°	0.56 mm (.022")	L/R	297-235RBX
In/Out Comp. .047"			

Package of 5



Full Patient Kit

READI-BASE®  
eXact™